

Science **made** smarter

Sera™ made mobile

The Sera™ Transport Cart offers great mobility and ease of use



Interacoustics

Audiometry

Tympanometry

ABR
www

OAE
www

Hearing Aid Fitting

Balance



Convenient transport

With the Sera™ Transport Cart, it is easier to transport Sera™ and its equipment, perfect for busy screening personnel. It stores the equipment in an organized way and lowers the risk of contamination by keeping the equipment away from the patient. Storing Sera™ on the Sera™ Transport Cart enables you to provide a hands-free service and extra care for your patients.

Features

- Five wheels for great mobility
- Handle that makes it easy to maneuver
- Retainers for all the equipment and accessories you might need
- Integrated storage bin for potential waste
- Storing bracket for Sera™ that secures it during transport while retaining ease of access for hand-held screening

Dimensions

- Weight including power supply: 11.3 kg
- Height from the floor to the top of the storage shelf: 1.16 m
- Wheelbase diameter: 0.53 m
- Utility bin: H 0.08 m × W 0.3 m × D 0.3 m
- Cart bin: H 0.08 m × W 0.36 m × D 0.23 m

Easy transportation of Sera™ and accessories



© Interacoustics a/s - 8520894 - 1 - 11/2019

